

VARROX® Evaporator

The Original Oxalic Acid Dihydrate Evaporator



Month of treatment

1	2	3	4	5	6	7	8	9	10	11	12
---	---	---	---	---	---	---	---	---	----	----	----

- Well tolerated by queen and bees
- High efficacy
- Easy dosage with measuring spoon

Best Application

- ✓ Brood-less colonies, November-February
- ✓ Ideal at temperature below 5 °C
- ✓ Use pure oxalic acid dihydrate for evaporation without residues.

Dosage

- ✓ 1 spoon makes 1 g
- ✓ For Langstroth (1 body) 1 g
- ✓ For Langstroth (2 bodies) and Dadant 2 g



1

1. Fill

Oxalic acid into the pan and place inside hive.



2

2. Connect

VARROX® to 12 V car battery.



3

3. Apply

inside tightly closed hive for 2 ½ minutes, disconnect from battery.



Andermatt
• • • • •
Canada

Andermatt Canada
Fredericton, N. B.
Tel. 506-444-5690
www.andermattcanada.com

0819

26.08.2019 09:18:18